

Dentosuc Spezial 2%

- Cleanser and disinfectant for dental suction systems-

Range of application:

Highly efficient disinfectant and cleanser for all kinds of dental suction systems. The regular use of Dentosuc Spezial 2% guarantees a hygienically clean suction system.

Effective range:

Bactericidal, fungicidal, virus inactivating (incl. HBV/HIV/HCV/BVDV)

Dosage and exposure time:

Acting time acc. DGHM/VAH:
2 % - 15 min.

Composition:

100 g solution contain
6.0 g dialkyldimethylammoniumchloride.

The tensides contained in the product are biodegradable acc. to regulation 648/2004 EC on detergents.

Environmental compatibility:

Dentosuc-Spezial 2% can be disposed of with the waste water as a diluted working solution after use. No remobilization of heavy metals, e.g. mercury, from amalgam due to a chlorine-free formulation.

CE
0124



Important notes:

Flammable. Flash point: acc. DIN 51755: 25°C. Application max. 50 ml per m². The total amount of application per room must not exceed 100 ml per m² of floor space. Keep container tightly closed. Do not breathe vapour/aerosol. Use only in well ventilated areas. No smoking. Keep away from sources of ignition! Keep out of reach for children

Marking acc. to Hazardous Substances Ordinance:

Irritating to eyes and skin. In case of contact with eyes, rinse immediately with plenty of water and seek medical advice. This material and its container must be disposed of in a safe way. Keep out of reach for children.



müller-
omicron
DENTAL

Müller-Omicron GmbH & Co. KG
Schlosserstr. 1, DE-51789 Lindlar
Tel. +49 22 66 4742-0

RFV 04/09